UD Smart®

Klinik

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Hibrit Protezler

Clinical

CASE-BOOK

Hybrid Prosthetics

UD UNITED DENTAL

İçindekiler | Contents

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5 UD Smart Nedir ?
What is UD Smart?

5 UD Smart Nedir ?
What is UD Smart?

8 Vaka 1:
Case 1:
Dr. Dt. Tolga Pekperdahcı &
Prof.Dr. İsmail Doruk Koçyiğit

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Dear Colleagues,

As we are all aware, the world has engaged in a huge transformation. Naturally, dentistry also has to take its share in this transformation.
This transformation has been changing not only profession's outlines but patients' expectations as well.
Such that, the patients, beyond health, are participating in aesthetic appearance, the material to be used and duration of the work to be done. Some these demands are appropriate. Whereas the other aspects should be considered by us professionals, but we don't have to power to resist the transformation.

At this point, we have to move with transformation and we have to shape the transformation together. As dentists, we responsible from aesthetic of prosthetic superstructures, right manufacturing, quality, and meeting patients expectations. Unfortunately, classic manufacturing techniques are insufficient in this transformation.
If "digital manufacturing techniques" that emerged as alternative are manufactured correctly, we can ensure development and transformation in the name of profession.
We have to confess that every digital manufacturing is not done correctly. Even worse, in some cases, classic manufacturing is better than digital manufacturing.

I congratulate United Dental that has brought all the truth together in digital manufacturing in our country. We have a lab that runs its own tests rather than believing in equipment suppliers' claims. This lab has done serious investments at huge costs to create a machine park to manufacture the highest quality. Apart from this, they have a patent of infrastructure which belongs to United Dental.

Patented product, UD Smart®, due to its superiorities is preferred by dentists in prothesis over implants.
Distributing impact evenly, passive fitting, extra lightness, and high resistance are the features make this patented technology superior.

As mentioned above, I congratulate United Dental for not manufacturing digital technology but manufacturing digital technology right. I hope that they will assist our profession for a long time..

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Nowadays with the advancing technology our options and opportunities, in the field of dental prosthesis have been increasing. Our patented product UD Smart®, has no better alternative for your cases which have hybrid prosthesis indication.

UD Smart® which has been brought about by the experience of our technicians and engineers is the fruit of many years of R&D.

In cases with minimum vertical dimension loss of 13,5 mm in single jaw in order to reduce the load of the implants and eliminate the complications such as broken screws, geometrical shapes like in the space frame technology specifically designed for every case are placed in the prosthesis, with this prosthesis gets lighter and stronger and at the same time with "hybrid manufacturing" method the "super fit" of the prosthesis to multiunit abutments is achieved.

The purpose of the casebook which was compiled from the cases of dentists from various regions of Türkiye to point out that it is possible for every dentist to finish such cases with the conventional and digital workflow which were perfected by us and to ensure that it plays, **UD Smart®** first in Türkiye then in the world scene.

In United Dental organization is our academic board which consists of internationally reputable dentists.

With our Academic Board, we provide the best service to our physicians both in the manufacturing and application of our products in the light of scientific developments in the sector.

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FP3, special to your cases, is a prothesis treatment method without alternative. Our expectation is to transmit vertical and horizon tal jaw movement impact which come onto prothesis on implants to each implant platform that are used as abutments in the same arch evenly. Load distribution can be done homogeneously and life span of implant and prothesis can be prolonged this way.

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
When compared with metal infrastructure prothesis obtained by lazer, milled from plate or by casting, UD SMART® which is superior to prothesis manufactured by these known methods is lighter in the ratio of 1⁄4. UD SMART®, even though manufactu- red from CoCr as all these known prothesis are, is the prothesis closest to the weight of natural teeth due to its inner structure.

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

If we are planning a screw restoration in your case, our expectation should be passive fit on implant abutment level. Whereas in super fit, there should not be a gap between two interpenetrating platforms as well as abutment's margin should meet abutment and prothesis' merging point. (The gap needed for cement restoration is not needed for screw restoration.) For this reason, the screw that binds abutment and prothesis will be centered and there will not be any impact on screw due to jaw movements. Therefore, due to super fit application, we will eliminate nightmares of screw cracks, neck region melting and neck region fractures. Super fit prothesis manufacturing can be done by more than one methods. The newest and the most precise method is Hybrid Manufacturing.

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UD SMART is manufactured by this method. Manufacturing the prothesis in layer manufacturing and milling the implant interface on the same prothesis in the CNC is called Hybrid Manufacturing. Namely, it is a method using two different manufacturing process in the same prothesis. Briefly, prothesis is being manufactu- red in 3D Metal Printer and interface on implant level or abutment level of the very same prothesis is milled in linear axed CNC up to 2 micron precision. In addition, the screw nest of the prothesis is milled with special T burs to allow dentist to obtain real torque in patient's mouth.

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
Image-3 : A vertical impact on full and motionless object (Osseointegrated implant) can be observed. As the level of impact increases, the object will be sprained and then cracked.

Image-4 : A horizontal impact on full and motionless object (Osseointegrated implant) can be observed. As the level of impact increases, the object will be sprained and then cracked.

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
The shape of the load does not carry importance, the rod pressures or pulls, distributing the entire load in the axle Therefore, compared to the full rods, the material has much more endurance in central and oblique loads

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Image-5 : Due to patented geometric structures in each mm2 of UD SMART®, the angle of impact, vertical or horizontal, does not matter any more. The stress caused by any impact on any part is pulled by other parts. Therefore, the load is distributed to the whole arch. UD SMART® has much more endurance than other inside-full-prothesis

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
*We provide possibilty to you dear dentists for a prothesis treatmet without alternative which could be used in All on 4 and All on X cases in various scenarios of patients with different physcial jaws, loss of extreme vertical dimension and bruxism. **UD SMART®** is designed specificaly to each patient.*

UD SMART® is super fit and manufactured with advanced engineering that can be covered up with porcelain and comes with a 7 year warranty. (images 6,7,8)

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 Our 53 year old male patient;
 There is no systemic disease.
 *In line with the **Anamnesis and Clinical Examination;***
 Current drug : None
 Alcohol consumption : 300 ml/week

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 Clinical examination revealed periodontal problems on each arch on almost all the teeth and a very obvious neglect of oral hygiene.
 In his clinical anamnesis, he has been treated several times due to his periodontal problems and he has his fixed porcelain crowns since 2000.
 He asked the treatment not to be edentulous in any stages of the treatment.
 *Initial situation of his upper jaw, hopeless teeth on both side of the **maxilla** and **mandible.***

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 Treatment Plan
 - *CBCT examination,*
 - *Extractions, upper and lower jaw teeth,*
 - *İmplant placement 15-12-22-25 and 45-42-32-35*
 - *Both Maxilla and Mandibula were immediately loaded with **Nucleoss T6 Implants,** multiunit abutments and temporary abutments*
 - *Four months later both upper and lower jaws treated with **UD Smart® Bar***
 - *Six months later post operative check*

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 Augmentin 1000 mg (2 x 1) was prescribed 24 hours prior to the surgery Surgery was proceeded in November 2020.
 *Local infiltration anesthesia with **Ultracain D-S Fort (Articain/epinefrine)** 5 X1.7 ml. Full thickness flap is elevated from the bone crest on the line of the attached gingiva. Teeth were extracted and all the surgical field cleaned from the **debritments, periodontal** wounds. Also **bone alveolplasty** was done with the burs and 2.5mm alveolectomy was performed due to the restorative space management.*

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 *4.5/14mm on the 15 position, 4.0/12 mm on the 12 position, 4.0/12 mm on the 22 position and 4.5/14mm on the 25th position were drilled and placed in the same time. Initial **torque values** was recorded and all of them were over **62**. On posterior sides **30° multiunit abutments** and on the anterior sides **0° multiunit abutments** were selected because of the immediate loading situation and fixed up to 15 Nm to the implants. After implants put into planned position, healing abutments are placed and flap closure is achieved with **4.0 PGA** (polyglicolic asid) suture.*

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 *After 15 days same surgical technic was used for lower jaw and implant placement. 4.5/12mm on the 45 position, 4.0/12 mm on the 42 position, 4.0/12 mm on the 32 position and 4.5/12mm on the 35th position were drilled and placed in the same time. Initial **torque values** was recorded and all of them were over **65**. On posterior sides 30° multiunit abutments and on the anterior sides 0° multiunit abutments were selected because of the immediate loading situation and fixed up to 15 Nm to the implants. After implants put into planned position, healing abutments are placed and flap closure is achieved with 4.0 PGA (polyglicolic asid) suture.*

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 *Preorganized prosthesis used as a **temporary prosthesis** with the proper technique. After United Dental lab procedures temporary prosthesis was fixed and occlusal check was finished. Soft diet and other postoperative instructions were given to the patient.*

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 *4 months later temporary prosthesis was removed, and clinical control was done. After uneventful healing **open tray impression** coping mounted on the multiunit abutments and impression was finished with the individual tray.*

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 *Stone model prepared and digitalized with the model scanner. After scanning procedure, the data was evaluated on the software and **UD Smart® Bar** to United Dental was designed and machined. In try-in session passive fit and porcelain process was completed by United Dental.*

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 Prosthesis screwed in and occlusal checks were completed, and routine check appointment was given to the patient.

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 Our patient states that the prosthesis raises her standard of living in her periodic controls, and she feels so comfortable and like her own teeth that she does not remember the presence of the prosthesis.

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
 Our 54 years old male patient;
 He does not have any systemic disease.
 *In line with the **Anamnesis and Clinical Examination;***
 Current drug : None
 Alcohol consumption : None
 Smoking habit : 1 pack a day

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 The upper jaw is completely edentulous, he has been using a full denture for about 2 years and the current prosthesis is incompatible. Two canine teeth in the mandible are periodontally hopeless and have caries. There is also a classical removable prosthesis made two years ago.
 She does not want to be toothless after the implant.

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 Treatment Plan
 - *CT examination*
 - *Extraction of lower canine teeth*
 - *Six implants were planned for regions 12-14-16-22-24-26, 42-44-46-32-34-36,*
 - *It is planned to install 6 **Nucleoss T6 Torq** implants in Maxilla and 6 **Nucleoss T6 implants** in Manbibula and load them immediately*
 - *After 4 months, both jaws will be rehabilitated with **UD Smart Bar***
 - *Post-operative control at 6 months*

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 *Augmentin 1000 mg (2 x 1) was prescribed 24 hours before surgery. Local infiltration anesthesia was performed with **Ultracain D-S Fort** (Articaine / epinephrine) 5 X 1.7 ml. The full thickness flap was opened from the crest to teeth 17-27. **Nucleoss T6 Torq** implants of 3.5/14mm, number 14 4.0/14mm, number 16 4.5/12mm, number 22 4.0/14mm, number 24 4.0/14mm and number 26 4.5/12 mm were placed in the tooth area 12. Initial torque values were recorded and all were 60 and above. Healing caps were placed after the implants were placed. The flap was closed with a 4.0 PGA (polyglycolic acid) suture. After the surgical procedure on the upper jaw, the lower jaw was started immediately.*

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 Nucleoss T6 implants of 3.5/12mm, size 44 4.1/12mm, number 46 4.1/12mm, number 32 3.5/12mm, number 34 4.1/12 and number 36 4.1/12mm were placed in the tooth area 42. Initial torque values were recorded and All of them were 65 and over. Healing caps were placed after the implants were placed. The flap was closed with a 4.0 PGA (polyglycolic acid) suture.

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 Impressions were taken with polyvinyl siloxane impression material using the open tray method from the upper and lower jaws immediately after the surgical procedures. After the models were obtained in the United Dental, the design in the exocad program was matched with the prostheses previously made to the patient, and the design in the form of the tooth that the patient liked was completed. This design was sent to a 3D dental printer and a resin temporary was produced. After laboratory procedures, a temporary prosthesis was inserted and occlusal control was completed. The patient was given a soft diet and other post-operative instructions.

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 *After 4 month, the temporary prostheses were removed and the necessary controls were made, and permanent measurements were taken with polyvinyl siloxane impression material using the method of connecting the implant level tray impression copings. After the models are prepared, they are monted on articulator with **Face-Bow** using temponary prostheses which the patient had been using.*

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 *Again, with reference to the temporary prostheses he used and said that he liked in general terms, after completion of the design on the exocad program the control prosthesis (Jig) produced with a **3D printer** was rehearsed in the patient's mouth.*

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 *On this rehearsal, the position, size and occlusion of the teeth were checked for the last time. After passive fit is checked and patient consent is obtained, it is sent to the United Dental Sent for **UD Smart® Bar** construction. The last changes were made and the bars were put into production.During the rehearsal, passive session and occlusion control were made and sent back to the United Dental for porcelain.*

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 After the treatment, we were informed that the patient was satisfied with the UD Smart hybrid prosthesis and that it was very comfortable to use.

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Our patient, who had implant losses in the past and was not satisfied with the roughness and weight of her old prosthesis, applied to our clinic, and implants were planned on the areas with sufficient bone on 3D tomography and applied intraoral.

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*Thanks to the digital workflow of UD dental, since the prosthetic connections of the previously made different brand implants and the newly made implants are different, thanks to the new generation software, while the **UD smart®** prosthesis is produced, the implant connections are also different, produced on the prosthesis.*

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Afterwards, the UD Smart® metal infrastructure, which was produced on the final prosthesis planning, was rehearsed both on the model and in the patient's mouth with a single screw test.

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After the dentin rehearsal, the final dentures of the patient were tightened with the recommended torques and delivered to the patient without any problems.

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
As a result, we completed our treatment with a lightweight, functional and more stable prosthesis by eliminating interconnections.

UD Smart®

Dünyanın En Hafif Protezi

The Lighest Prothesis Of The World

Doğala En Yakın Hibrit Protez
The closest hybrid prothesis to the natural

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 Our 60 years old female patient;
 He does not have any systemic disease.
 *In line with the **Anamnesis and Clinical Examination;***
 Current drug : None
 Alcohol consumption : None
 Smoking habit : None

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 The UD Smart® metal infrastructure, which was produced on the final prosthesis to United Dental planning, was rehearsed both on the model and in the patient's mouth with a single screw test.

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 Atrophic jaws were rehabilitated with zygomatic implant and mandibular All on 4 implant hybrid prostheses.

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*Our 60 year old male patient applied to our clinic with pain in his teeth, swelling on his face and bleeding gums. In the clinical and radiographic examinations, deep periodontal pockets and periodontitis and periodontal abscesses characterized by severe horizontal and angular bone losses were detected in the existing teeth. **Implant-abutment** and **abutment-prosthetic** restoration around the implants made to the patient at different times **Peri-implatitis**, which is caused by maladaptation in their relationship and characterized by bone loss, was detected.*
Although it is possible to save some of the teeth in the mouth from a periodontal perspective, it is not possible to make a fixed prosthetic restoration that the patient desires in this way. Therefore, it was decided to remove all teeth and old implants in the mouth and to perform an over-implant prosthetic restoration with new implants.

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*All the teeth and implants in the mouth were removed in a single session under local anesthesia, the hard tissue around them was cleared and the tissues were curetted. finely developed **periodontitis** and **peri-implantitis** extensive bone defects, graft from the **"creos™ xenogain"** type and **"creos™" xenoprotect”** was augmented with the collagen membrane and membrane fixing pins.After 8 months of waiting for new bone growth to occur, 5 implants were placed in the upper jaw and 4 implants in the lower jaw **(Nobel Biocare Parallel™ CC)**. A torque higher than 45 Ncm was obtained for all implants. Thus, in the same session, the **“Multi-unit abutment plus”** and **“Temporary Snap Coping Multi-unit Plus”** abutments and the lower and upper total prostheses with acrylic teeth, which we had previously prepared and checked for occlusion, were attached to the implants to be used as fixed temporary prosthesis. At the end of 4 months, measurements were taken with the impression technique with an open troy.*

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 *In order for the prosthesis to be prepared to sit passively on all implants in the patient's mouth, firstly, controls were made in the patient's mouth with **"plaster jigs"**, which are an indicator of the accuracy of the measurements we took. When there was no break in the **"plaster jigs"**, metal infrastructure preparations were started. It is necessary to increase the vertical dimension of the prosthesis in order to regain the lost vertical dimension, especially in prostheses made in cases of total edentulism, which develops as a result of advanced bone destruction. This causes an increase in the weight of the prosthesis.*

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 *Although titanium can be used as a substructure for lightweight prostheses, there are clinical problems with the connection of **titanium-porcelain**. For this reason, **UD Smart® (United Dental Dental Prosthesis Laboratory Şişli-Istanbul)** prosthesis production technique, which is produced with the Hybrid Manufacturing method, is used for both the reliability of working on porcelain on CoCr alloy and lightness. Thus, metal supported porcelain prosthesis was produced by United Dental by combining both the advantages of laser sinter production technique and the precision of **CNC** engraving systems.*

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*Metal supported porcelain prostheses produced with By United Dental the **UD Smart®** technique were screwed to the "multi unit abutment plus" in the patient's mouth with a torque of 15Ncm and the screw holes in the prosthesis were closed first with teflon tape and then with composite restorations.*

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 Our 65 year old female patient;No general systemic problems. Clinical and tomographic examination; The patient had tooth loss in both the lower and upper jaws, chronic lesions around the teeth, and advanced bone loss and mobility in the existing teeth due to chronic periodontitis.

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 After tooth extraction under surgery and local anesthesia, cleaning of the lesions, 4 implants (Nobel Biocare Parallel CC) with 2 standard and 2 angles in the upper jaw with the concept of All-on-4, 4 flat implants in the lower jaw (Nobel Biocare Parallel CC) and 1 standard total of 5 implants were placed, including straight implants (AlphaBio Neo). By choosing appropriate multi units, immediate prostheses were made on both jaws on the same day with acrylic temporary prostheses with the help of multi unit temporary titanium abutments.

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 *Permanent Prosthesis; the patient was able to come to start his permanent prosthesis after about 8 months due to his job. A pattern resin jig was prepared by **United Dental** using open impression posts on the model obtained from the impression taken during surgery. Cut jigs were tied back in the mouth and open measurements were taken with prepared troys. The next day, the plaster jig prepared for the impression control was checked in the mouth and it was determined that there was no breakage and the impression was uneventful. In the next step, a rehearsal was performed with resin rehearsal prostheses prepared with the guidance of the patient and the physician. In this rehearsal, the final prosthesis was tested in terms of both aesthetics, lip support and occlusion. After obtaining the patient's consent, necessary occlusal adjustments were made. Afterwards, the **UD Smart®** metal infrastructure, which was produced on the final prosthesis planning, was rehearsed both on the model and in the patient's mouth with a single screw test.*

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The patient has been using his prostheses for 2 years without any problems, very happy with both aesthetics and function.

Dr. Dt. Emre Çimen
Çene Cerrahisi ve İmplantoloji
Maxillofacial surgery and Implantology

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 Our 58 year old male patient; No general systemic problems. In accordance with the clinical and tomographic examination, tooth loss in both the lower and upper jaws, chronic lesions around the teeth, and advanced bone loss and mobility in the existing teeth due to chronic periodontitis were observed.

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 After tooth extraction and cleaning of the lesions under surgery and local anesthesia, 2 standard and 2 angled 4 implants (Nobel Biocare Parallel CC) were placed in the upper jaw with the All on 4 concept, and 6 standard implants (AlphaBio Neo) in the lower jaw. By choosing appropriate multi units, immediate prostheses were made on both jaws on the same day with United Dental acrylic temporary prostheses with the help of multi unit temporary titanium abutments.

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 *Permanent Prosthesis; After a 4 month waiting period, a pattern resin jig was prepared by **United Dental** using open impression posts on the model obtained from the impressions taken during surgery. Cut jigs were tied back in the mouth and open impression were taken with prepared troys. The next day, the plaster jig prepared for the measurement control was checked in the mouth and it was determined that there was no breakage and the impression was uneventful. In the next step, a rehearsal was performed with resin rehearsal prostheses prepared with the guidance of the patient and the dentist. In this rehearsal, the final prosthesis was tested in terms of both aesthetics, lip support and occlusion. After obtaining the patient's consent, necessary occlusal adjustments were made. Afterwards, the **UD Smart®** metal infrastructure, which was produced on the final prosthesis planning, was rehearsed both on the model and in the patient's mouth with a single screw test.*

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
About 6 months after the delivery of the prosthesis, chipping occurred on the porcelain surface due to the patient's teeth grinding habit. With the advantage of making a screw retained, the prosthesis was removed and sent to United Dental, where it was quickly repaired, and a night guard was applied to the patient. Except for this patient-related event, our patient has been using his prostheses for about 2 years without any problems, being quite happy with both aesthetics and function.

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 Age and gender of the patient : 65 years / male
 General systemic details: Regular drug intake – No
 Smoking : No

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 A brief description of the case: There are broken teeth in the lower and upper jaw region, teeth with periapical lesions and teeth with periodontal problems. There was alveolar bone loss.
 Clinical indication: Extraction of all upper and lower teeth and guided bone regeneration.

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 *12 Astra Ev implants were applied to our case, 6 in the lower jaw and 6 in the upper jaw. In the upper jaw of this case, we worked with the occlusal screw hybrid manufacturing interface system and the CoCr interface engnoving system. In the lower jaw we applied the **UD Smart®** system. Screw retained was performed in both upper and lower prostheses and we ensured passive fit. After checking the accuracy of the manual impressions with a plaster jig, we provided the smile design and occlusion control by using resin milling. Afterwards, we completed our case by moving to the production stages of our passively seated screw retained screw prostheses. We checked the fit of our prostheses with X-ray.*

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 *Thoughts on **UD Smart®**:*
 Being lighter is a big advantage. It can be preferred especially in coses where vertical dimension is high. Although the vertical distance in the lower jaw was greater, the lower jaw became a lighter prosthesis.

Dr. Dt. Emre Çimen
Çene Cerrahisi ve İmplantoloji
Maxillofacial surgery and Implantology

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 Our 46 year old male patient; No general systemic problems. Clinical and tomographic examination; tooth loss in both the lower and upper jaw of the patient, chronic lesions around the teeth, and chronic lesions in the existing teeth.There was advanced bone loss and mobility due to periodontitis.

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 After tooth extraction and cleaning of the lesions under surgery and local anesthesia, 2 standard and 2 angled 4 implants (Nobel Biocare Parallel CC) were placed in the upper jaw with the All-on-4 concept, and 6 straight implants (AlphaBio Neo) in the lower jaw. By choosing appropriate multi units, immediate prostheses were made on both jaws on the same day with acrylic temporary prostheses with United Dental the help of multi unit temporary titanium abutments.

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*A pattern resin jig was prepared by **United Dental** using open impression posts on the model obtained from the impression taken during surgery. Cut jigs were tied back in the mouth and open measurements were taken with prepared troys. In the next step, a rehearsal was performed with resin rehearsal prostheses prepared with the guidance of the patient and the physician. After obtaining the patient's consent, necessary occlusal adjustments were made. Afterwards, the **UD Smart®** metal infrastructure, which was produced on the final prosthesis planning, was rehearsed both on the model and in the patient's mouth with a single screw test.*

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After the dentin rehearsal, the final dentures of the patient were tightened with the recommended torques and delivered to the patient without any problems. The patient has been using his prostheses for 1 year without any problems, being very happy with both aesthetics and function.

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
A 45-year-old male patient applied to our clinic to get information about his treatment because he did not like his existing teeth aesthetically. In the radiography and clinical examination, bone loss was observed in the posterior region due to long-term edentulism, along with the loss of periodontal support tissue. After the patient was informed that it was not possible to perform a fixed treatment with his existing teeth, the treatment options were explained to our patient (after our patient did not want sinus lift/bone augmentation surgeries in terms of time), it was decided to continue with the "All on 4" treatment option.

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After the extraction of the existing teeth in the mouth, implant loading (DXL, Gelsenkirchen, Germany) was performed immediately. After the primary stability was taken, fixed immediate loading from the lower maxilla total prosthesis was performed on the 2nd day by using the "multi unit temporary abutments" of the company DXL. Our patient was informed about the care of the prostheses that he will use for 3 months.

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*After the completion of the osseointegration (average 3 months), the permanent prosthesis was made. An open impression was taken using the multiunit open troy impression post. Controlling the passive fit in all screwed prostheses is important for the long-term success of the implant and prosthesis. After the determination of the vertical dimension, after the passive fit control with the plaster jig, the infrastructure rehearsal stage was started. Due to the placement of the maxillary angled implants close to the canine region of our patient, a prosthesis was planned for the upper jaw with a narrow occlusal table and the amount of cantilever was determined according to the AP distance in **EXOCAD** during the design phase in order to ensure that the final prosthesis was finished with a short arc concept. The metal interior in the upper jaw was designed with space frame technology and **UD Smart®** Smart Hybrid Prosthesis produced with Hybrid Manufacturing technique and screw retained metal supported prosthesis in the lower jaw was preferred.*

Msc. Dt. Gökhan Tunalıgil
Ağız Diş ve Çene Cerrahisi, Oral İmplantoloji
Oral and Maxillofacial Surgery, Oral Implantology

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 Total extraction and upper All on 4 and lower conventional implant treatment were planned in a 52 year old male patient due to periapical infections and bone loss.

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 *In order to make the upper All on 4 prosthesis lighter, it was decided to make it with the **UD Smart®** system.*

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 At the end of the 3 month osseointegration period, the patient's fixed prostheses were applied.

Uzm. Dr. Dt. Kerem Turgut Atasoy
Ağız Diş ve Çene Cerrahisi
Oral and Maxillofacial Surgery

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In a 60 year old female patient who had old implants in the lower jaw and could not reach the prosthetic parts, implant removal and All on 4 and upper tooth extraction and All on X implant treatment were planned.

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*At the end of the 3-month waiting period in the patient without immediate loading prepaped by United Dental **UD Smart®** hybrid prostheses were applied.*

Uzm. Dr. Dt. Kerem Turgut Atasoy
Ağız Diş ve Çene Cerrahisi
Oral and Maxillofacial Surgery

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Total extraction and fixed prosthetic treatment with zygoma implants upper jow and standard dental implants ower jow were planned in a 56 year old female patient due to severe periodontal destruction and bone tissue deficiencies.

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Intraoral view of the patient before the impression

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The position, size and closure of the teeth were checked with a resin proof prepared with a 3D printer in accordance with the closure we received on the Exocad program at United Dental

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*During the three-month recovery period, the prosthesis was not loaded on the patient, and at the end of three months, fixed prostheses were applied to the patient with the **Ud Smart** system and delivered after the impression taken with the open troy impression method.*

Uzm. Dr. Dt. Kerem Turgut Atasoy
Ağız Diş ve Çene Cerrahisi
Oral and Maxillofacial Surgery

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 *63 year old male patient. All on 4 in the upper jaw, conventional dental implant treatment in the lower jaw and **UD Smart®** hybrid fixed prostheses were planned due to partial missing teeth and pneumatized maxillary sinuses.*

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
 *Measurements were taken with an open spoon measuring system at the 3rd month from the implants that were not immediately loaded. In the 3D printer trynins, the patient's criteria such as occlusion, smile line, and midline were evaluated. Produced successfuly at United Dental produced **UD Smart®** hybrid prostheses were attached to multi unit abutments with passive fit.*

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 In a 77 year old retired nurse patient, removal of dental implants made in an external center, tooth extraction and fixed hybrid prosthesis supported by lower and upper implants were planned.

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
 Dental implants were placed without flaps in the upper jaw and open surgery in the lower jaw, and the 3-month recovery period was overcome without wearing a temporary prosthesis.

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
 *At the end of the treatment, **UD Smart®** prostheses were fitted.*

Msc. Dt. Sanlı Volkan Öz
Oral İmplantoloji ve Çene Cerrahisi
Oral Implantology and Maxillofacial Surgery

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
 Our 45 year old female patient;
 There is no systemic disease.
 *In line with the **Anamnesis and Clinical Examination;***
 Current drug: None
 Alcohol consumption: None

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 *In the clinical and radiographic examinations, deep periodontal pockets and periodontitis and periodontal abscesses characterized by severe horizontal and angular bone losses were detected in the existing teeth. **Implant-abutment** and **abutment-prosthetic** restoration around the implants made to the patient at different times **Periimplatitis**, which is caused by maladaptation in their relationship and characterized by bone loss, was detected.*

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 Although it is possible to save some of the teeth in the mouth from a periodontal perspective, it is not possible to make a fixed prosthetic restoration that the patient desires in this way. Therefore, it was decided to remove all teeth and old implants in the mouth and to perform an over-implant prosthetic restoration with new implants.

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Metal supported porcelain prostheses produced with the UD Smart® technique by United Dental were screwed to the "multi unit abutment plus" in the patient's mouth with a torque of 15Ncm and the screw holes in the prosthesis were closed first with teflon tape and then with composite restorations.

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 Our 58 year old female patient;
 There is no systemic disease.
 In line with the ***Anamnesis and Clinical Examination;***
 Current drug : None
 Alcohol consumption : None

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 At the end of the 3-month osseointegration period, fixed prostheses were applied to the patient at United Dental.

Uzm. Dr. Dt. Kerem Turgut Atasoy
Ağız Diş ve Çene Cerrahisi
Oral and Maxillofacial Surgery

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 60 year old male patient had old ball-attachment dental implants that had bone loss in the upper jaw and had to be removed, and teeth that had to be extracted due to bone loss in the lower jaw.

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 *While zygoma implant was applied in the left upper region of the patient, conventional implant treatment was planned in the lower jaw and **UD Smart®** hybrid prosthesis was planned.*

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 *Fixed prosthetic treatment was completed with **UD Smart®** hybrid prostheses applied at the end of the three-month recovery period.*

Uzm. Dr. Dt. Kerem Turgut Atasoy
Ağız Diş ve Çene Cerrahisi
Oral and Maxillofacial Surgery

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 In a 60 yea old patient residing abroad, All on 4 dental implant treatment was planned by performing tooth extractions and zygoma implants on the upper left and conventional dental implants on the lower.

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
 *During the 3-month healing period, the implants were not loaded and at the end of 3 months, the patient's **UD Smart** ®hybrid prostheses were applied.*

Msc. Dt. Gökhan Tunalıgil
Ağız Diş ve Çene Cerrahisi, Oral İmplantoji
Oral and Maxillofacial Surgery, Oral Implantology

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 In a 55 year old male patient, vertical bone loss and dehiscence in the buccal tissue were observed in previous dental implants in the mandible.

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 It was decided to remove the old dental implants and apply dental implants to the lower jaw again.

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 *At the end of the 3-month recovery period, the patient's impression were taken with an open tray impression system and passive fit was confirmed with plaster jigs. **UD Smart®** hybrid prosthesis was attached to the lower jaw of the patient.*

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We have been serving many European conutries such as United Kingdom Germany, France and Belgium.

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You are very close to obtaining fast and effective results via special customer representative who is specialized in his/her field.

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
Digital impression delivery options:
send your files via Wetransfer.com Use info@udlabor.com as recipient email.

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UD/ST® UNITED DENTAL SMART TECHNOLOGY

UD/R&D® UNITED DENTAL RESEARCH & DEVELOPMENT

UD/DW® UNITED DENTAL DIGITAL WORKFLOW

UD/TS® UNITED DENTAL TRAINING SERVICE

UD® UNITED DENTAL Academy

www.uniteddental.com.tr

UD20230120-NS